


NOVOSOUND


Wearables Evaluation Kit



**laptop not included*


Novosound Wearables Evaluation Kit: a hardware and software platform for measuring and monitoring of soft tissue (skin, muscle, vascular, etc).

Kit includes a **flexible wearable sensor** comprising a central high amplitude transmit element and dual 16 element broad bandwidth receive arrays.



This array is designed to be connected to Novosound's **wireless wearable ultrasound platform**.

Proprietary electronics include a single channel transmit and 32 channel receive pulser system, with WiFi connectivity to stream data to a customer's own display unit*



www.novosound.net

info@novosound.net